

HumaCLIA 150

HumaCLIA SR Immunoassay Menu



- Ready-to-use reagents for immediate use, no waiting time for temperature equilibration
- Small kits with 2 x 50 tests per cartridge, ideal also for labs with smaller sample numbers
- Calibrators in the reagent kits, for streamlined logistics and improved workflow efficiency in the lab
- Fully-validated applications, ensuring reliable result quality

Thyroid	REF
TSH	84600
fT3	84620
fT4	84610
T3	84615
T4	84605
anti-TG	84625
anti-TPO	84630
TG*	84635

Tumor Marker	REF
PSA	83130
fPSA	83125
AFP	83100
CEA	83105
CA-125	83110
CA 19-9	83120
CA 15-3	83115

Fertility hormones	REF
AMH	82150
HCG	82100
FSH	82120
LH	82115
Progesterone	82105
Testosterone	82110
Prolactin	82125
Estradiol (E2)	82130
DHEA-S*	82145

Cardiac Marker	REF
CK-MB	85040
Troponin I	85030
BNP*	85045
Myoglobin	85035

Bone Marker	REF
PTH**	85010
25-OH Vitamin D	85000

Anemia	REF
Ferritin**	85049
Folate*	85060
Vitamin B12	85055

Inflammation	REF
PCT	85020

Metabolism	REF
Cortisol*	85070

Infectious diseases	REF
anti-TP (Syphilis)*	81100

* Coming soon

** Non-CE – only for use outside Europe

HumaCLIA 150

Controls and Consumables



HumaCLIA SR Controls	REF	To be used with HumaCLIA SR
Immunoassay Multi Control <i>2 levels each with 2 x 3 ml</i>	84850	TSH, ft3, ft4, T3, T4, anti-TG, anti-TPO, TG, AFP, CEA, CA-125, CA 19-9, CA 15-3, FSH, LH, HCG, Prolactin, Progesterone, Testosterone, Estradiol (E2), DHEA-S, PSA**, fPSA**, 25-OH Vitamin D, PTH**, Ferritin**, Folate, Vitamin B12, Cortisol, Troponin I, CK-MB, Myoglobin
AMH Control <i>2 levels each with 2 x 3 ml</i>	82850	AMH
PCT Control <i>2 levels each with 2 x 2 ml</i>	85820	PCT
Anti-TP Control* <i>2 levels each with 2 x 2 ml</i>	81800	anti-TP (Syphilis)

HumaCLIA SR Diluent	REF	To be used with HumaCLIA SR
Sample Diluent B <i>2 x 45 ml</i>	80110	HCG

HumaCLIA SR Consumables	REF
Reaction Vessels <i>2 x 1,000 pcs</i>	15910/11
Pre-Trigger Solution <i>4 x 220 ml</i>	80050
Trigger Solution <i>4 x 220 ml</i>	80055
Wash Buffer <i>4 x concentrate for 10 l each</i>	80060



Scan the QR-code or visit
www.human.de/cliareagents
for latest product availability.

* Coming soon

** Non-CE – only for use outside Europe